

## **Product Data Sheet**

## Bowe Systec/Protechno Card Attachment Labels

PRODUCT SERIES	1) BOW — 21
	2) BOW — 21R - WITH PRINTED REGISTRATION MARK
	3) PRO - 14
MATERIAL — CARRIER	1) 76.2 μ (.003") CLEAR POLYOLEFIN
	2) 76.2 μ (.003") CLEAR POLYOLEFIN
	3) 127.0 μ (.005") WHITE COMPOSITE
MATERIAL - LINER	1) 81.0 μ (.0032") DIFFERENTIAL RELEASE, WHITE GLASSINE
	2) 81.0 μ (.0032") DIFFERENTIAL RELEASE, WHITE GLASSINE
	3) 112.0 μ (.0044") DIFFERENTIAL RELEASE, POLY LAMINATED WHITE KRAFT
SIZE	17mm X 38mm (0.66" X 1.5")
CHARACTERICTIC	A LIVERIA LARGE MATERIAL THAT RELEADED FACILIVE ROMELIE LINER CONCIDENT
CHARACTERISTIC	A HYBRID LABEL MATERIAL THAT RELEASES EASILY FROM THE LINER. CONSISTENT ADHESION TO PVC. POLYESTER, POLYCARBONATE AND MULTI COMPOUNDED CARD
	SUBSTRATES. LEAVES NO ADHESIVE RESIDUE UPON REMOVAL FROM CARD.
MINIMUM ORDER QUANTITY	- 10 ROLLS
SHELF LIFE/STORAGE	ALL DIAGENICS CARD SERVICES INC. LABEL PRODUCTS MAY BE SAFELY STORED FOR A PERIOD OF SIX (6) MONTHS AND MUST BE STORED IN ITS PROTECTIVE COVER, IF SUPPLIED. LABELS MUST BE KEPT IN A CLEAN AREA, FREE FROM EXCESSIVE MOISTURE & SHOULD BE STORED AT 40 - 60% RELATIVE HUMIDITY, AWAY FROM ANY DIRECT HEAT SOURCE OR SUNLIGHT. AMBIENT TEMPERATURE SHOULD BE BELOW 38°C.

The information contained relates only to the specific indentified substrate materials. This information should be used as a guide only in determining the product suitability for the application, and should be confirmed by independent test procedures before use.